

Atys medical Basic 3.2 et 3.4



Caractéristiques	BASIC L	BASIC 1	BASIC 2	BASIC 3.1	BASIC 3.2	BASIC 3.4
Doppler 4 & 8 MHz		✓	✓	✓	✓	✓
Remplissage veineux bilatéral	✓		✓	✓	✓	✓
PPG bilatérale	✓		✓	✓	✓	✓
PVR bilatérale					✓	✓
Pression systolique Doppler				✓	✓	✓
Pression systolique PPG					2	4
IPS Doppler				✓	✓	✓
IPS PPG					✓	✓
IPSO					✓	✓
Test de marche				✓	✓	✓
Vitesse onde de pouls					✓	✓

PPG= Photo-pléthysmographie
PVR = Pneumo-pléthysmographie