

ECHOGRAPHE EDAN U-60 Couleurs



Sonde convexe C352UB - Frequencies: 2.5, 3.5, 4.5, H2.5, H2.7 MHz
Applications: Abdomen, OB, GYN, et Urology
Angle: 90°



Sonde linéaire L1042UB - Frequencies: 6.5, 7.5, 8.5, H4.5, H4.7 MHz
Applications: Small Parts, Vascular, Musculoskeletal, Superficia
Largeur 40 mm



Sonde micro-convexe C422UB - Frequencies: 5.5, 6.5, 7.5, H4.5, H4.7 MHz
Applications: Cardiac and Pediatric.
Angle: 90°



Sonde endo-vaginale E612UB - Frequencies: 5.5, 6.5, 7.5, H4.5, H4.7 MHz



Sonde linéaire L552UB - Frequencies: 6.5, 7.5, 8.5, H4.5, H4.7 MHz
Applications: Musculoskeletal, Vascular, Breast, and Orthopedics.
Scanning Length: 60 mm.



Imprimante Sony UP-897MD

L'imprimante Sony UP-897MD est l'imprimante la plus petite et la plus légère de sa catégorie. Cette imprimante compacte est l'imprimante la plus rapide pour imprimer des impressions en format A6 en noir et blanc. Le résultat des impressions est nette et de haute qualité.



Mobile Trolley pour EDAN

